



OCTG Series,

OCT-LKx, and OCT-RAx

**OCT Standard Scanner,
Scan Lens Kit,
and Standard Scanner
Reference Arm Adapter**

User Manual



















Original User Manual – not translated

Table of Contents

Chapter 1	Warning Symbol Definitions.....	1
Chapter 2	Introduction	2
2.1.	Safety	2
2.2.	Care and Maintenance	3
2.2.1.	Optical Cleaning	3
2.2.2.	Service	3
2.2.3.	Accessories and Customization	3
Chapter 3	Scanner Compatibility.....	4
Chapter 4	Installation	5
4.1.	OCTG Mounting	5
4.2.	OCTG Connections	6
4.2.1.	Connecting the Electrical Control Interface.....	6
4.2.2.	Connecting the Optical Fiber.....	7
4.3.	Integration	8
Chapter 5	Description	9
5.1.	Theory	9
5.1.1.	Signal Generation	9
5.1.2.	Limitations	9
5.2.	Optical Design	10
5.2.1.	Common Path Setup	10
5.2.2.	Dual Path Setup (OCTG-NR).....	10
5.2.3.	Realization.....	11
5.3.	Components	12
5.3.1.	OCTG Base Module	13
5.3.2.	OCTG Reference Module	15
5.3.3.	OCT Scan Lens Kit (Accessory)	16
5.3.4.	Reference Arm Adapter (Accessory)	17
5.3.5.	Dove Tail Mount	18
5.4.	Dimensions	19
Chapter 6	Troubleshooting.....	20
Chapter 7	Certifications and Compliance.....	21
Chapter 8	Specifications.....	22
8.1.	Reflectivity Scanning Mirror	23
Chapter 9	Warranty.....	24
9.1.	Imaging Systems	24
9.2.	Non-Warranty Repairs	24
9.3.	Warranty Exclusions	24
Chapter 10	Regulatory	25
10.1.	Waste Treatment is Your Own Responsibility	25
10.2.	Ecological Background	25
Chapter 11	Thorlabs OCT Support Contact.....	26

Chapter 1 Warning Symbol Definitions

Below is a list of warning symbols you may encounter in this manual or on your device.

Symbol	Description
	Direct Current
	Alternating Current
	Both Direct and Alternating Current
	Earth Ground Terminal
	Protective Conductor Terminal
	Frame or Chassis Terminal
	Equipotentiality
	On (Supply)
	Off (Supply)
	In Position of a Bi-Stable Push Control
	Out Position of a Bi-Stable Push Control
	Caution: Risk of Electric Shock
	Caution: Hot Surface
	Caution: Risk of Danger
	Warning: Laser Radiation
	Caution: Spinning Blades May Cause Harm

Chapter 2 Introduction



2.1. Safety

Please read this manual carefully before operating the OCTG standard scanner. Please also read any manuals for the systems being connected to the OCTG.

All statements regarding safety and technical specifications will only apply when the unit is operated correctly.

 **ATTENTION** 

This equipment is intended for laboratory use only and is not certified for medical applications, including, but not limited to, life support situations.



 **WARRANTY WARNING** 

There are sensitive electronic and optical parts in the OCTG.

Any modification or servicing of this system by unqualified personnel renders Thorlabs free of any liability.

Any modification of the galvanometer scanners or the camera may cause loss of the factory optical alignment.

This device can only be returned for service when it is packed into the complete original packaging, including all foam packing inserts. Please contact Thorlabs' OCT support (see Chapter 11) for replacement packaging if the original packaging has been lost.

 **LASER RADIATION WARNING** 

When a light source (e.g. SLD, laser) is being coupled into the OCTG, please observe the appropriate laser safety precautions for your own protection. The appropriate laser safety precautions depend on the light source coupled into the OCTG.

During normal operations, laser light will be present within the scanner and will also be emitted from the scanner. Also laser light may be emitted from unexpected locations, such as if the fiber has been disconnected from the body or if the reference arm has been disconnected.

In addition, the OCTG is an optical system that can influence the divergence of the beam. This can cause a change of laser class of the light source, especially if the OCTG is used without an objective.

Always turn off the light source before changing or adjusting the OCTG configuration or accessories as the objective, lens kit, reference arm adapter, or sample z-spacer. For Thorlabs OCT base units, turn off the OCT base unit main power to turn off the light source.

2.2. Care and Maintenance

The system should be treated with care, particularly during transportation and unpacking. Hitting or dropping the system can damage the unit and lower system performance. If mishandling occurs, misalignment of the optical components may occur, leading to a decrease in image quality. In this situation, the system should be realigned by qualified personnel. Do not store or operate in a damp, closed environment.

- Do not store or operate on surfaces that are susceptible to vibrations.
- Do not expose to direct sunlight.
- Do not use solvents on or near the equipment.
- Keep the unit away from dust, dirt, and airborne contaminants, such as cigarette smoke. The system is not designed for outdoor use. Protect the equipment from rain, snow, and humidity.
- Do not subject the equipment to mechanical and thermal extremes. Protect the equipment from rapid variations in temperature.
- Handle all electrical and fiber connectors with care. Use of excessive force to form electrical or fiber connections may damage the connectors.

2.2.1. Optical Cleaning

The most common cause of low signal intensity is dirtying of the fiber due to airborne contaminants. To minimize the fiber's exposure to air, avoid unnecessary disconnections of the optical fiber patch cable. Ensure that the connection is tight, and keep the fiber as straight as possible without placing it under tension. It is also advisable to check the fiber when making other adjustments to the optical system, such as changing the objective.

Thorlabs' Fiber Inspection Scope (Item # FS200) can help determine when the fiber needs cleaning. We recommend our Fiber Connector Cleaner (Item # FCC-7020) for quickly cleaning the fiber tips.

2.2.2. Service

Only trained and approved Thorlabs personnel are allowed to service the system. Please contact Thorlabs' OCT support (see Chapter 11) for more information.

2.2.3. Accessories and Customization

The OCTG series standard scanners are Thorlabs-qualified accessories for Thorlabs' OCT Systems (i.e., our CALLISTO, GANYMEDE, VEGA and TELESTO). We strongly suggest using Thorlabs' OCT-LKx scan lens kits as well as OCT-RA reference arm adapter kits with the OCTG as they were specifically designed to work together.

In order to achieve the intended performance, this scanner should only be used with qualified parts. Please hold a conversation with Thorlabs' OCT support (see Chapter 11) to determine if other parts you wish to use are compatible. Any modification or servicing of this system by unqualified personnel renders the warranty null and void, leaving Thorlabs free of any liability.

Chapter 3 Scanner Compatibility

The OCTG standard scanner is a standalone, preassembled, integrated scanner intended for the use together with a Thorlabs OCT base unit such as the CALLISTO, VEGA, GANYMEDE or TELESTO series.

This scanner is available in two different setups

- The OCTG-xxxNR standard scanner is designed for dual path setups using a dedicated external reference to create the interferometric signal. These scanners does NOT comprise an interferometer.
- The other OCTG-xxx scanner is designed for common path setups and comprise an interferometer consisting of a beam splitter and a reference arm creating the interferometric signal. The OCTG-xxx is available in two versions for different wavelength ranges.
 - The OCTG-900 for OCT systems working in the 900 nm regime.
 - The OCTG-1300 for OCT systems working in the 1300 nm regime.

In this manual we will use abbreviation for the OCTG scanner as follows:

- OCTG OCTG-900; or OCTG-1300

For common statements the abbreviation “OCTG” is used for both setups.

The OCTG scanners are fully compatible with all Thorlabs OCT base units of the CALLISTO, VEGA, GANYMEDE, and TELESTO series.

The table below gives a short overview of the different standard scanners, their usable wavelength range and lists preferred OCT base units.

Standard Scanner	Wavelength Range	OCT Base Unit
OCTG-900	850 nm – 1000 nm	CALxxx GANxxx
OCTG-1300	1200 nm – 1400 nm	TELxxx
OCTG-1300NR	1200 nm – 1400 nm	VEGxxx

Table 1 Usable Wavelength Range of OCT Scanner

The most selective optical component is a mirror mounted on the Y galvanometric scanner of the scanner set.

A reflectivity graph of the scanning mirror is given in Figure 21.

Detailed information about the spectral performance of the different mirrors are available upon request. Please contact Thorlabs' OCT support (see Chapter 11) for details

Chapter 4 Installation

The OCTG standard scanner is a standalone, preassembled, integrated accessory to an OCT base unit. It should be securely mounted to an optical table or breadboard with minimal vibrations. We recommend mounting the OCTG to a Thorlabs OCT-Stand.

4.1. OCTG Mounting

To mount the OCTG in the OCT-Stand, gently slide the dovetail of the OCTG into the slide of the OCT-Stand.




Figure 1 Mounting the OCTG in the OCT-Stand

4.2. OCTG Connections

4.2.1. Connecting the Electrical Control Interface

Attach the included electrical control cable to the OCTG. You may use either side of the cable since the plugs are identical. The OCT scanner's electrical control interface is located at the top of the OCTG base module. Align the red dot of the plug to the alignment mark of the port.

Push the connector into the receptacle until a "click" sound is heard. This click indicates that the connector is locked.




Figure 2 Installation of the Control Cable Connector at the OCTG

The remaining plug of the electrical control cable will be attached to the OCT base unit. Make sure that the base unit is switched off.

The connection is located at the rear of the base unit. For installation, align the red dot, facing the alignment mark in the base unit. Push the connector into the receptacle marked "Control" until a "click" sound is heard. This click indicates that the connector is locked.

4.2.2. Connecting the Optical Fiber



ATTENTION



When installing the fiber, make sure that the fiber tip does not get contaminated by dust. Thorlabs' Fiber Inspection Scope (Item # FS200) and Fiber Connector Cleaner (Item # FBC1) are useful for keeping the optical path clean. Do not touch the fiber tip!

Attach the optical fiber to the OCTG, as illustrated in Figure 3 below. Either end of the fiber patch cable may be used to connect to the OCTG. Remove the dust cap from each fiber end and store them with the system packaging. The FC/APC fiber connection is located at the top of the OCTG base module side by side with the electric interface. Insert the fiber tip into the center bore of the fiber connection, then secure the tip by gently rotating the locking cap clockwise.

The fiber connector needs to be oriented such that alignment key slides into the key slot of the fiber connector (as shown in Figure 3 below). If the key is NOT properly aligned with respect to the key slot, you will still be able to screw in the fiber connector, but significant light intensity loss and focal shift will result from this incorrect connection.




Figure 3 Installation of the Fiber at the OCTG Standard Scanner

4.3. Integration

For the full integration into such a system, please refer to the user manual of the OCT base unit.



Figure 4 Fully integrated Callisto System with 900 nm Standard Scanner (Item # OCTG-900) Mounted on a Stand (Item # OCT-Stand) with Translation Stage (Item # OCT-XYR1) operated by a Base Unit (Item # CAL110)

Chapter 5 Description

5.1. Theory

5.1.1. Signal Generation

Spectral domain optical coherence tomography (SD-OCT) generates cross-sectional images up to several millimeters deep into tissue. The images are assembled by performing a series of scans at adjacent, increasing depths, allowing 2D and 3D reconstruction of the specimen.

Therefore the interference between the light coming back from the specimen and from a reference is sampled for a broad range of different wavelengths. This is performed either by using a broadband light source divided into separate wavelength using a spectrometer or by a laser with small bandwidth quickly tuned over a broad wavelength range.

The phase delay of the back-reflected and back-scattered light (with respect to the stationary reference) is recorded as a function of wavenumber, and a Fast Fourier Transform (FFT) yields the cross-sectional images as a function of sample depth.

5.1.2. Limitations

The spatial resolution and sensitivity of the OCTG scanner depends on several parameters, including the following:

- **Wavelength Range:** The optical components within the OCTG are optimized for a specific wavelength range, depending upon the model. For the usable wavelength range please refer to Table 1.
- **Optical Power:** The sensitivity of the OCTG is directly related to the intensity of the light returning from the sample. Factors that can reduce the collected light intensity from the sample fiber include: dirty fibers, blocked/cropped beams, and condensed water in the environment.
- **Physical Movements:** OCT systems use a camera to detect the phase relation of the light returning from the sample. Even small movements of the specimen in relation to the optical reference arm can "wash out" the wavenumber-resolved phase contrast, affecting the image.
- **Imaging:** In a fiber-based OCT setup, the light returning from the sample is focused into the core of an optical fiber. Hence, the fiber can be thought of as a spatial filter for the light. This filter has an effective diameter, referred to as the "mode field diameter. For single mode propagation, mode field diameter is larger than the fiber core diameter. Poor focusing, caused by optical aberration or misalignment, therefore leads to loss of contrast and sensitivity.

5.2. Optical Design

5.2.1. Common Path Setup

As shown below in Figure 1, in this configuration, the scanner is factory-configured such that the sample beam and reference beam are generated after the beam leaves the fiber. This allows you to use single-mode optical fiber to transport the beam into the scanner while minimizing the use of free-space propagation.

This approach avoids problems that can degrade image quality, related to optical phenomena like polarization mode dispersion (PMD) and birefringence, and makes the performance of the system independent from the length of the single-mode fiber.




Figure 5 Diagram of the Common Path Setup

5.2.2. Dual Path Setup (OCTG-NR)

Shown below in Figure 6 is the beam geometry when the scanner is ordered in a special configuration without reference. The sample beam and reference beam are generated within different fibers, before the beam exits into free space. In this configuration, the scanner becomes the sample arm of the interferometer. By using two different fibers, the beamsplitter cube used in the Common Path Setup is no longer needed.

This configuration allows single mode optical fiber to be used to a greater extent within the setup. While this approach is able to provide greater sensitivity due to the absence of the beamsplitter cube (which reduces the intensity of the light that returns to the fiber), it is significantly more sensitive to the optical phenomena mentioned before. Please contact Thorlabs' OCT support (see Chapter 11) for details.




Figure 6 Diagram of the Dual Path Setup

5.2.3. Realization

The basic optical layout of the OCTG scanner in common path layout is illustrated below in Figure 7.




Figure 7 Optical Layout OCTG

The output of an FC/APC fiber is collimated and routed to a beam splitter cube. Here the beam is divided into a sample beam and a reference beam, similar to a Michelson interferometer. The sample beam is routed over two galvanometer actuated mirrors to allow for scanning in two axes. The scan objective then focuses the beam in the sample. Back-scattered and back-reflected light is collected by the scan objective again and travels back to the fiber. The light reflected into the reference arm is retro-reflected back into the fiber. There is an optimum intensity for the reference light that can be adjusted using the reference intensity adjustment knob which will open or close the variable aperture inside the OCTG.

In the dual path configuration of the OCTG-NR, the reference path components are not included.




Figure 8 Optical Layout OCTG Common Path

5.3. Components

The OCTG Standard Scanner are to be used in combination with Thorlabs' OCT base units.



Figure 9 OCT Modules OCTG and OCTG-NR

The major components of the OCTG Standard Scanner are an OCTG base module and a Reference Arm Module if applicable. We recommend also considering other accessories like the OCT-LKx scan lens kit and the OCT-RAx reference arm adapter.

The following sections describe the two modules and possible accessories in detail.

ATTENTION

Always turn off the light source before changing or adjusting OCTG configuration or accessories such as the objective, scan lens kit, reference arm adapter, or sample z-spacer. For Thorlabs OCT base units, turn off the OCT base unit main power to turn off the light source.

Please contact a member of the Thorlabs' OCT support team to determine if other parts you wish to use are compatible (see Chapter 11). Any modification or servicing of this system by unqualified personnel renders the warranty null and void, leaving Thorlabs free of any liability.

5.3.1. OCTG Base Module

The OCTG base module provides high-speed, two-dimensional (X and Y) raster scans of the specimen. The clear aperture of the scan mirrors used within is \varnothing 6 mm. The module also contains a high-resolution video camera for recording the sample during the measurement.




Figure 10 OCT Base Module

The electrical control interface hosts the included Thorlabs electric control cable to connect the scanner with a Thorlabs OCT base unit. Please contact Thorlabs' OCT support (see Chapter 11) for information regarding the pin configuration.

The optical fiber port is a FC/APC receptacle.

The optical output to the reference arm is equipped with a circular iris aperture which is manipulated using the reference adjustment knob. In order to adjust the reference intensity, it may be necessary to rotate the reference length adjustment knob. Rotating clock-wise increases the reference intensity, rotating counter clock-wise decreases it.



ATTENTION



The iris aperture might be damaged by inadequate torque. In positions “open” and “closed” an enlarged mechanical resistance indicates the limit of the travel range.

In order to adjust the reference intensity adjustment knob, pull the knob approximately 5 mm outwards until you feel the knob coming to a rest.



Figure 11 Reference Intensity Adjustment

As a qualitative indication, observe the reference intensity bar in the OCT software. Please refer to the Software Manual for additional guidance.

5.3.2. OCTG Reference Module

The reference module contains a mounted mirror (i.e., retro reflector) that reflects the beam from the light source back into the OCT interferometer.



Figure 12 OCTG Reference Module

In order to match the optical path length in this reference arm to the optical path length of the light from the sample, it may be necessary to translate the mirror along the axis of the optical system. Length adjustments can be performed by rotating the reference length adjustment knob.



Figure 13 Reference Length Adjustment

Rotating clock-wise increases the reference path length, rotating counter clockwise reduces it.

The position of the mirror can be monitored using the reference length indicator on the front side of the module.

The full adjustment range of the reference module is >12 mm while the intended standard position is around 3 mm giving an adjustment range of -2 mm / +10 mm.

5.3.3. OCT Scan Lens Kit (Accessory)

The OCT scan lens kit from Thorlabs are specially designed accessories to support telecentric scanning over a wide range and to compensate the dispersion mismatch of the scanner. The OCT-LKx are fully compatible with the OCTG.

If you ordered a scan lens kit and reference adapter with the OCTG, these items are pre-installed.

If not already installed the installation procedure follows these three steps.

- Insert the included OCT scan lens by screwing it in,
- Inserting the illumination module
- Secure the illumination module using the four included M2.5x6 cap screws




Figure 14 Installation of Scan Objective and Illumination

- Another part of the lens kit is the dispersion compensation set. This has to be inserted into the reference adapter OCT-RAx



Figure 15 Dispersion Compensation Set



ATTENTION



**Optics are sensitive components!
When handling the objective lens take care to avoid touching or harming the optical surfaces.**

Changing the objective of the OCTG requires changing the probe configuration data in the software package ThorImage OCT as well.

This is performed automatically using the calibration procedure built-in the ThorImage OCT Software package version 4.1.4 and higher. A detailed description of this procedure is given in the ThorImage OCT Operating Manual.

5.3.4. Reference Arm Adapter (Accessory)

To match the required optical path length and the dispersion for the used OCT scan lens kit the installation of the reference arm adapter OCT-RAx is strongly suggested.

If you ordered a reference arm adapter with the OCTG, this item is pre-installed.

If not already installed the dispersion compensation kit provided with the lens kit and must be installed

- Insert the dispersion compensation kit, which is part of the lens kit, into the reference adapter
- Fix it using the SM05 retaining ring




Figure 16 Installation of Dispersion Compensation Set

For the installation procedure of the OCT-RAx follow these steps.

- Remove the OCT reference module by loosening the three fixed M3x6 cap screws
- Insert the reference arm adapter onto the OCT base module
- Inserting the OCT reference module onto the reference arm adapter
- Secure the assembly using the three included M3 cap screws



Figure 17 Installation of Reference Arm Adapter

5.3.5. Dove Tail Mount

The OCTG ships with a dovetail mount at the back side (as seen in Figure 18 below). This accessory allow the scanner to be mounted in standard Thorlabs OCT-Stand.



Figure 18 OCTG-1300 with Dove Tail Mount

The OCT-Stand is a dedicated stand for OCTG scanners.

The OCTG scanners are attached to the focus block of the OCT-Stand using a spring loaded mount accepting the dove tail mount of the OCTG. The focus block can be rotated 360° around the Ø1.5" rod, and features 30 mm of travel with fine and coarse adjustment knobs.




Figure 19 OCT Stand with OCT-XYR1

5.4. Dimensions

The dimensions of the OCTG series standard scanners are given in the following drawing. All dimensions are in mm.




Figure 20 Drawing showing Dimensions, Mounting Options, and Focal Plane

The opto-mechanical specifications of the OCT scan lens kits are listed in the table below.

Scan Lens Kit #	No OCT-LK	OCT-LK2 (-BB)	OCT-LK3 (-BB)	OCT-LK4 (-BB)
Field of View	-	6 mm x 6 mm	10 mm x 10 mm	16 mm x 16 mm
Mechanical Height (H) OCTG	139 mm / 71.5 mm	142.5 mm / 71.5 mm	162 mm / 71.5 mm	189 mm / 71.5 mm
Barrel Height (B)	-	25.5 mm	24.0 mm	37.0 mm
Working Distance (WD)	-	3.4 mm	24.9 mm	41.6 mm

Table 2 Data for OCT Scan Lens Kits

Changing the objective of the OCTG requires changing the probe configuration data in the software package ThorImage OCT as well.


This is performed automatically using the calibration procedure built-in the ThorImage OCT software package version 4.1.4 and higher. A detailed description of this procedure is given in the ThorImage OCT operating manual.

Chapter 6 Troubleshooting

Problem	Possible Cause	Recommended Solution
Poor Reference Light Intensity	Fiber Not Connected	Remove and Reconnect Fiber, Ensuring that Alignment Key is Inserted into Key Slot
	Aperture is Too Small	Open Aperture
	Fiber Tip is Dirty	Clean Fiber Tip (Thorlabs' MCC-7020 Fiber Connector Cleaner Recommended)
	Other Reason	Contact OCT Service (See Chapter 11)
No Image is Obtained	Optical Path Length of Reference and Sample Arms is Not Matched	Adjust Reference Arm Length
	Beam is Blocked	Check for Obstructions in Optical Path
	USB Cable is Loose	Reconnect USB Cable
	PC Crashed	Restart PC
	Other Reason	Contact OCT Service (See Chapter 11)
Low Scan Resolution	Dispersion in Reference and Sample Arms is Not Matched	Check if Dispersion compensation is installed
Bad Image Quality	Image Obtained is Being Mirrored	Adjust the Distance Between the Objective and the Sample. The Image Should Move Towards the Top of the Computer Window
	Optical Path Length of Sample Arm is Too Short	Move Sample Away From Objective
	Reference Intensity is Too High or Too Low	Close or Open Aperture Iris to Adjust Intensity
	Other Reason	Contact OCT Service (See Chapter 11)
Flipped Image	Optical Path Length of Reference Arm is Incorrect	Adjust Reference Arm Length

Table 3 Troubleshooting

Chapter 7 Certifications and Compliance



THORLABS
www.thorlabs.com

EU Declaration of Conformity

in accordance with EN ISO 17050-1:2010

We: **Thorlabs GmbH**
Of: **Hans-Boeckler-Str. 6, 85221 Dachau/München, Deutschland**

in accordance with the following Directive(s):

2006/42/EC	Machinery Directive (MD)
2014/30/EU	Electromagnetic Compatibility (EMC) Directive
2011/65/EU	Restriction of Use of Certain Hazardous Substances (RoHS)

hereby declare that:
Model: **OCTG-xxxx(NR)**

Equipment: **OCTG Standard Scanner**


is in conformity with the applicable requirements of the following documents:

EN ISO 12100	Safety of Machinery. General Principles for Design. Risk Assessment and Risk Reduction	2010
EN 61326-1	Electrical Equipment for Measurement, Control and Laboratory Use - EMC Requirements	2013

and which, issued under the sole responsibility of Thorlabs, is in conformity with Directive 2011/65/EU of the European Parliament and of the Council of 8th June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, for the reason stated below:


does not contain substances in excess of the maximum concentration values tolerated by weight in homogenous materials as listed in Annex II of the Directive

I hereby declare that the equipment named has been designed to comply with the relevant sections of the above referenced specifications, and complies with all applicable Essential Requirements of the Directives.

Signed:  On: **15 February 2017**

Name: **Dorothee Jennrich**
Position: **General Manager**

EDC - OCTG-xxxx(NR) -2017-02-15



Chapter 8 Specifications

OCTG	
Optical Specifications	
Center Wavelength	900 nm, or 1300 nm
Clear Aperture	Ø6 mm (Max)
Reference Length Fine Adjustment	-2 mm / +10 mm
Scan Distance (Objective Shoulder)	15.0 mm / 28 mm
General Specifications	
Video Camera	Color CMOS
Weight of Scanner	1 kg (2.2 lbs)
Storage / Operating Temperature	10 °C to 35 °C
Dimensions of OCTG (L x W x H)	85 mm x 60 mm x 139 mm
Dimensions of OCTG-NR (L x W x H)	85 mm x 60 mm x 71 mm
Airborne Noise Emission	< 70 dBA

Table 4 Specifications OCTG

OCT Scan Lens Kits			
Objective Item #	OCT-LK2(-BB)	OCT-LK3(-BB)	OCT-LK4(-BB)
Field of View	6 mm x 6 mm	10 mm x 10 mm	16 mm x 16 mm
Barrel Height (B)	25.5 mm	24.0 mm	37.0 mm
Working Distance	3.4 mm	24.9 mm	41.6 mm

Table 5 Specifications Scan Lens Kits

8.1. Reflectivity Scanning Mirror

The OCTG series standard scanners are equipped with a semitransparent galvo mirror to enable video camera imaging.*




Figure 21 Reflectivity of Scanning Mirror

* In systems delivered prior to February 2017 a coating with different characteristics was used. For further information about their performance, please contact Thorlabs' OCT support (see Chapter 11).

Chapter 9 Warranty

9.1. Imaging Systems

Thorlabs offers a one-year warranty on the OCTG standard scanner.

9.2. Non-Warranty Repairs

Products returned for repair that are not covered under warranty will incur a standard repair charge in addition to all shipping expenses. This repair charge will be quoted to the customer before the work is performed.


9.3. Warranty Exclusions

The stated warranty does not apply to products which are (a) specials, modifications, or customized items (including custom patch cables) meeting the specifications you provide; (b) ESD sensitive items whose static protection packaging has been opened; (c) items repaired, modified, or altered by any party other than Thorlabs; (d) items used in conjunction with equipment not provided by or acknowledged as compatible by Thorlabs; (e) subjected to unusual physical, thermal, or electrical stress; (f) damaged due to improper installation, misuse, abuse, or storage; (g) damaged due to accident or negligence in use, storage, transportation, or handling.

Chapter 10 Regulatory

As required by the WEEE (Waste Electrical and Electronic Equipment Directive) of the European Community and the corresponding national laws, Thorlabs offers all end users in the EC the possibility to return “end of life” units without incurring disposal charges.

- This offer is valid for Thorlabs electrical and electronic equipment:
- Sold after August 13, 2005
- Marked correspondingly with the crossed out “wheelie bin” logo (see right)
- Sold to a company or institute within the EC
- Currently owned by a company or institute within the EC
- Still complete, not disassembled and not contaminated



Wheelie Bin Logo

As the WEEE directive applies to self-contained operational electrical and electronic products, this end of life take back service does not refer to other Thorlabs products, such as:

- Pure OEM products, that means assemblies to be built into a unit by the user (e. g. OEM laser driver cards)
- Components
- Mechanics and optics
- Left over parts of units disassembled by the user (PCB's, housings etc.).

If you wish to return a Thorlabs unit for waste recovery, please contact Thorlabs or your nearest dealer for further information.

10.1. Waste Treatment is Your Own Responsibility

If you do not return an “end of life” unit to Thorlabs, you must hand it to a company specialized in waste recovery. Do not dispose of the unit in a litter bin or at a public waste disposal site.

10.2. Ecological Background

It is well known that WEEE pollutes the environment by releasing toxic products during decomposition. The aim of the European RoHS directive is to reduce the content of toxic substances in electronic products in the future.

The intent of the WEEE directive is to enforce the recycling of WEEE. A controlled recycling of end of life products will thereby avoid negative impacts on the environment.

Chapter 11 Thorlabs OCT Support Contact

If you have a technical question or issue on Thorlabs OCT products, please refer directly to the OCT Support team located in Luebeck, Germany.

OCT Support

Thorlabs GmbH
Maria-Goeppert-Straße 9
23562 Lübeck
Germany
Tel: +49-(0)8131-5956-0
Fax: +49-(0)8131-5956-99
www.thorlabs.de
Email: oct-support@thorlabs.com

Chapter 12 Thorlabs Worldwide Contacts

For technical support or sales inquiries, please visit us at www.thorlabs.com/contact for our most up-to-date contact information.



USA, Canada, and South America

Thorlabs, Inc.
sales@thorlabs.com
techsupport@thorlabs.com

Europe

Thorlabs GmbH
europe@thorlabs.com

France

Thorlabs SAS
sales.fr@thorlabs.com

Japan

Thorlabs Japan, Inc.
sales@thorlabs.jp

UK and Ireland

Thorlabs Ltd.
sales.uk@thorlabs.com
techsupport.uk@thorlabs.com

Scandinavia


Thorlabs Sweden AB
scandinavia@thorlabs.com

Brazil

Thorlabs Vendas de Fotônicos Ltda.
brasil@thorlabs.com

China

Thorlabs China
chinasales@thorlabs.com



THORLABS
www.thorlabs.com