

## H101-MCL-Z100/500

**Stage compatibility: fits on Mad City Labs Z100/Z500 piezo Z-stages**

**The chamber can fit:**

- 35/60 mm Petri
- 1x3 in. / 1x2 in. chambered slides
- 6-12-24-48-96 MW plates

**Features:**

- Uniformly heated by means water circulation in water tight channels present in the main body and in the lid
- Magnets allow easy interchange chamber insert
- Magnetic locks hold Petri, Slides in correct position inside the chamber
- Perfusion holes: 12 channels for 2.5 mm O.D. tubings are available in chamber riser
- Removable riser to fit standard multiwell plates or to perform perfusion