

H101-NIKON-NZ100/200/500-N

Stage compatibility: fits on Ti-S-E or Ti-S-ER + Mad City Labs Piezo Z-stages

The chamber can fit:

- 35/60 mm Petri
- 1x3 in. / 1x2 in. chambered slides
- 6-12-24-48-96 MW plates

Features:

- Uniformly heated by means water circulation in water tight channels present in the main body and in the lid
- Magnets allow easy interchange chamber insert
- Magnetic locks hold Petri, Slides in correct position inside the chamber
- Perfusion holes: 12 channels for 2.5 mm O.D. tubings are available in chamber riser
- Removable riser to fit standard multiwell plates or to perform perfusion