

H101-NIKON-TI-SR

Stage compatibility: fits on Nikon XY stages: Nikon T/Ti-SR

The chamber can fit:

- 35/60 mm Petri
- 1x3 in. / 1x2 in. chambered slides

Features:

- Uniformly heated by means water circulation in water tight channels present in the main body and in the lid
- Magnets allow easy interchange chamber insert
- Magnetic locks hold Petri, Slides in correct position inside the chamber
- Perfusion holes: 8 channels for 2.5 mm O.D. tubings are available in chamber riser