

## H101-OLYMPUS-IX-SUSP

**Stage compatibility: fits on Olympus XY stages: IX-SUSP, GX-SVR, GX-SFR, BX3-SSU, IX2-SFR, IX2-SVL-2**

**The chamber can fit:**

- 35/60 mm Petri
- 1x3 in. / 1x2 in. chambered slides

**Features:**

- Uniformly heated by means water circulation in water tight channels present in the main body and in the lid
- Magnets allow easy interchange chamber insert
- Magnetic locks hold Petri, Slides in correct position inside the chamber
- Perfusion holes: 8 channels for 2.5 mm O.D. tubings are available in chamber riser