

H201-PRIOR-NZ100-H107

Stage compatibility: fits on Prior NZ100/200 in Prior H107

The chamber can fit:

- 35/60 mm Petri
- 1x3 in. chambered slides
- No insert is required

Features:


- The chamber hosts #1 35mm petri or #1 chamber slide
- Flexible chamber height thanks to the removable chamber riser required to perform perfusion
- Perfusion holes: 12 channels for 2.5mm O.D. tubings are available in chamber riser








H201-PRIOR-NZ100-H107

H201-PRIOR-NZ100-H107 - 3D View



H201-PRIOR-NZ100-H107

H201-PRIOR-NZ100-H107 - Components



H201-PRIOR-NZ100-H107 - Dimensions