





Ø 1.38in
[35.0mm]





M32 x 0.75
EXTERNAL THREAD

M25 x 0.75
INTERNAL THREAD
(24.18mm CLEAR APERTURE)



FOR INFORMATION ONLY
NOT FOR MANUFACTURING PURPOSES

DRAWING PROJECTION					THORLABS www.thorlabs.com	
	NAME	DATE		M32x0.75 EXTERNAL TO M25x0.75 INTERNAL THREAD ADAPTER		
DRAWN	DW	05/AUG/16		MATERIAL BRASS		
APPROVAL	PN	05/AUG/16				
COPYRIGHT © 2016 BY THORLABS				REV A		
VALUES IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUNDOFF ERRORS				ITEM # M32M25S		APPROX WEIGHT .02 kg