




D2N FEMALE DOVETAIL




IMAGING PLANE



Ø 1.48in [37.6mm]

D1N MALE DOVETAIL



FOR INFORMATION ONLY  
NOT FOR MANUFACTURING PURPOSES

DRAWING  
PROJECTION



	NAME	DATE
DRAWN	TB	25/MAY/16
APPROVAL	PN	25/MAY/16

COPYRIGHT © 2016 BY THORLABS

VALUES IN PARENTHESIS ARE CALCULATED  
AND MAY CONTAIN ROUND OFF ERRORS

**THORLABS**  
www.thorlabs.com

TRINOCULARS WITH EYEPIECE & IR FILTER

MATERIAL

N/A

REV  
A

ITEM #

WFA4000

APPROX WEIGHT

1.72 kg