



FOR INFORMATION ONLY  
NOT FOR MANUFACTURING PURPOSES

DRAWING PROJECTION			<b>THORLABS</b> www.thorlabs.com	
	NAME	DATE	Fluorescent Calibration	
DRAWN	GY	21/JAN/19	Microscope Slide, Green, 1" x 3"	
APPROVAL	PN	21/JAN/19	MATERIAL	REV
COPYRIGHT © 2019 BY THORLABS			Cast Acrylic	A
VALUES IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS			ITEM #	APPROX WEIGHT
			FSK2	.002 kg