



THIS SET OF MICROSCOPE SLIDES INCLUDES 5 RED SLIDES
 ALL SLIDES HAVE THE SAME DIMENSIONS
 PLASTIC SLIDE BOX INCLUDED



FOR INFORMATION ONLY
 NOT FOR MANUFACTURING PURPOSES

DRAWING PROJECTION			THORLABS www.thorlabs.com	
	NAME	DATE	Fluorescent Calibration Microscope Slides, Red, Pack of 5	
DRAWN	GY	21/JAN/19	REV A	
APPROVAL	PN	21/JAN/19	MATERIAL	Cast Acrylic
COPYRIGHT © 2019 BY THORLABS			ITEM #	APPROX WEIGHT
VALUES IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS			FSK6-P5	0.06 kg