



SOLD IN PACKS OF 1000

NOTES AND SPECIFICATION:

1. THICKNESS TOLERANCE: $\pm 0.005\text{mm}$.
2. DIMENSIONAL TOLERANCE: $\pm 0.2\text{mm}$.
3. REFRACTIVE INDEX: 1.523 AT 589.3nm.
4. AVERAGE TRANSMISSION: (d=0.15mm) $91.7\% \pm 0.1\%$.

FOR INFORMATION ONLY
NOT FOR MANUFACTURING PURPOSES

DRAWING PROJECTION			THORLABS www.thorlabs.com	
	NAME	DATE	COVER GLASS #1.5, 22x22mm	
DRAWN	GY	16/JUL/19	MATERIAL	
APPROVAL	PN	16/JUL/19	SCHOTT D263	REV A
COPYRIGHT © 2019 BY THORLABS			ITEM #	APPROX WEIGHT
VALUES IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS			CG15CH	0.21g