



NOTES AND SPECIFICATION:

1. THICKNESS TOLERANCE:  $\pm 0.005\text{mm}$
2. DIMENSIONAL TOLERANCE:  $\pm 0.2\text{mm}$
3. REFRACTIVE INDEX: 1.523 AT 589.3nm
4. AVERAGE TRANSMISSION: (d=0.15mm) 91.7%  $\pm 0.1\%$

**SOLD IN PACKS OF 100**

FOR INFORMATION ONLY  
NOT FOR MANUFACTURING PURPOSES

DRAWING PROJECTION			<b>THORLABS</b> www.thorlabs.com	
	NAME	DATE	Cover Glass #1.5H, 24x50mm	
DRAWN	GY	17/JUL/19	MATERIAL	
APPROVAL	PN	17/JUL/19	SCHOTT D263 M	REV A
COPYRIGHT © 2019 BY THORLABS			ITEM #	APPROX WEIGHT
VALUES IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS			CG15KH1	0.51g