



NOTES AND SPECIFICATION:

1. THICKNESS TOLERANCE: $\pm 0.005\text{mm}$
2. DIMENSIONAL TOLERANCE: $\pm 0.5\text{mm}$
3. REFRACTIVE INDEX: 1.523 AT 589.3nm
4. AVERAGE TRANSMISSION: (d=0.15mm) $91.7\% \pm 0.1\%$

SOLD IN PACKS OF 100

FOR INFORMATION ONLY
NOT FOR MANUFACTURING PURPOSES

DRAWING PROJECTION			THORLABS www.thorlabs.com	
	NAME	DATE	COVER GLASS #1.5H, Ø 12mm	
DRAWN	GY	17/JUL/19	MATERIAL	
APPROVAL	PN	17/JUL/19	SCHOTT D263 M	REV A
COPYRIGHT © 2019 BY THORLABS			ITEM #	APPROX WEIGHT
VALUES IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS			CG15NH1	0.05g