



SOLD IN PACKS OF 1000

NOTES AND SPECIFICATION:

1. THICKNESS TOLERANCE: $\pm 0.005\text{mm}$.
2. DIMENSIONAL TOLERANCE: $\pm 0.5\text{mm}$.
3. REFRACTIVE INDEX: 1.523 AT 589.3nm.
4. AVERAGE TRANSMISSION: (d=0.15mm) $91.7\% \pm 0.1\%$.

**FOR INFORMATION ONLY
NOT FOR MANUFACTURING PURPOSES**

DRAWING PROJECTION			 www.thorlabs.com	
	NAME	DATE		
DRAWN	GY	17/JUL/19	COVER GLASS #1.5, 25mm DIA	
APPROVAL	PN	17/JUL/19	MATERIAL	REV
COPYRIGHT © 2019 BY THORLABS <small>VALUES IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS</small>			SCHOTT D263	A
			ITEM #	APPROX WEIGHT
			CG15XH	0.21g