



NOTES AND SPECIFICATION:

1. THICKNESS TOLERANCE: $\pm 0.005\text{mm}$
2. DIMENSIONAL TOLERANCE: $\pm 0.2\text{mm}$
3. REFRACTIVE INDEX: 1.523 AT 589.3nm
4. AVERAGE TRANSMISSION: (d=0.15mm) $91.7\% \pm 0.1\%$

SOLD IN PACKS OF 1000

FOR INFORMATION ONLY
NOT FOR MANUFACTURING PURPOSES

DRAWING PROJECTION			THORLABS www.thorlabs.com	
	NAME	DATE	Cover Glass #0, 22x22mm	
DRAWN	GY	17/JUL/19		
APPROVAL	PN	17/JUL/19	MATERIAL	REV
COPYRIGHT © 2019 BY THORLABS			SCHOTT D263 M	A
VALUES IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS			ITEM #	APPROX WEIGHT
			CG00C	0.2g