



NOTES AND SPECIFICATION:

1. REFRACTIVE INDEX: 1.523 AT 589.3nm
2. AVERAGE TRANSMISSION (d=1mm) 91.2%±1%
3. COLOR OF MARKING AREA: WHITE
4. 90° GROUND EDGES

SOLD IN PACKS OF 100

FOR INFORMATION ONLY
NOT FOR MANUFACTURING PURPOSES

DRAWING PROJECTION			THORLABS www.thorlabs.com	
	NAME	DATE	76x26mm slide, 1mm thick, positively charged	
DRAWN	AG	7/AUG/17		
APPROVAL	PN	11/AUG/17	MATERIAL	REV
COPYRIGHT © 2017 BY THORLABS			SODA LIME GLASS	A
VALUES IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS			ITEM #	APPROX WEIGHT
			MS10PC1	4.65g