



SOLD IN PACKS OF 1000

NOTES AND SPECIFICATION:

1. REFRACTIVE INDEX: 1.52 AT 589.3nm.
2. AVERAGE TRANSMISSION (d=1mm) 91.2% ± 1%.
3. COLOR OF MARKING AREA: WHITE.
4. 90° GROUND EDGES.

FOR INFORMATION ONLY
NOT FOR MANUFACTURING PURPOSES

DRAWING PROJECTION			THORLABS www.thorlabs.com	
	NAME	DATE	MICROSCOPE SLIDE, 1mm THICK UNIMARK WHITE	
DRAWN	AS	24/NOV/15		
APPROVAL	PN	25/NOV/15	MATERIAL	REV
COPYRIGHT © 2015 BY THORLABS			SODA LIME GLASS	A
VALUES IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS			ITEM #	APPROX WEIGHT
			MS10UW	4.65g