



SOLD IN PACKS OF 200

NOTES AND SPECIFICATION:

1. REFRACTIVE INDEX: 1.523 AT 589.3nm
2. AVERAGE TRANSMISSION (d=1mm) 91.2% ± 1%
3. COLOR OF MARKING AREA: WHITE
4. 90° GROUND EDGES

FOR INFORMATION ONLY
NOT FOR MANUFACTURING PURPOSES

DRAWING PROJECTION			THORLABS www.thorlabs.com	
	NAME	DATE	MICROSCOPE SLIDE, 1mm THICK UNIMARK WHITE	
DRAWN	AG	7/AUG/17	MATERIAL SODA LIME GLASS	REV A
APPROVAL	PN	11/AUG/17		
COPYRIGHT © 2017 BY THORLABS			ITEM #	APPROX WEIGHT
VALUES IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS			MS10UW2	4.65g