



- NOTES AND SPECIFICATION:
1. REFRACTIVE INDEX: 1.52 AT 589.3nm.
  2. AVERAGE TRANSMISSION (d=1mm) 91.2%±1%.
  3. CAVITY: 1 CAVITY Ø 15-18mm, DEPTH 0.6-0.8mm.
  4. 90° GROUND EDGES.

**SOLD IN PACKS OF 50**

FOR INFORMATION ONLY  
NOT FOR MANUFACTURING PURPOSES

DRAWING PROJECTION			<b>THORLABS</b> www.thorlabs.com	
NAME	DATE	MICROSCOPE SLIDE, 1.35mm THICK, ONE CAVITY		
DRAWN GY	4/AUG/2020	MATERIAL	SODA LIME GLASS	
APPROVAL PN	4/AUG/2020	ITEM #	MS15C1	APPROX WEIGHT 6.18g
COPYRIGHT © 2020 BY THORLABS		VALUES IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS		
		REV B		