



**SOLD IN PACKS OF 50**

**NOTES AND SPECIFICATION:**

1. REFRACTIVE INDEX: 1.52 AT 589.3nm.
2. AVERAGE TRANSMISSION (d=1mm) 91.2% ± 1%.
3. CAVITY: 2 CAVITY Ø 15-18mm, DEPTH 0.6-0.8mm.
4. 90° GROUND EDGES.

**FOR INFORMATION ONLY  
NOT FOR MANUFACTURING PURPOSES**

DRAWING PROJECTION			 <a href="http://www.thorlabs.com">www.thorlabs.com</a>	
	NAME	DATE		
DRAWN	AS	24/NOV/15	MICROSCOPE SLIDE, 1.5mm THICK TWO CAVITIES	
APPROVAL	PN	25/NOV/15		
COPYRIGHT © 2015 BY THORLABS			MATERIAL	REV
VALUES IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS			SODA LIME GLASS	A
ITEM #		APPROX WEIGHT		
MS15C2		6.20g		